

Certificate of Authenticity

Shaded Cercon® Zirconia

The dental restorations used for your case were constructed of FDA approved 510k materials and hand crafted in an ISO 13485:2003 certified dental laboratory.

This is your assurance of quality and safety in materials and construction.

Item	Product Name	Manufacturer	Country of Origin	FDA 510k	Elements & Compounds
Zirconium oxide core	Cercon	Dentsply	Germany	K013230	ZrO ₂ 93%, Y ₂ O ₃ 5%, Hf ₂ O ₃ <2%, other oxides <1%
Coloring Liquid	ZirkonZahn Ice	ZirkonZahn	Germany	K061851	Colored Salts: PEG, FeCl ₃ , MnCl ₂ , CrCl ₂
Porcelain	Vintage ZR	Shofu Dental	USA	K052033	Ceramic powder made from SiO ₂ , Al ₂ O ₃ , K ₂ O, Na ₂ O, CaO, B ₂ O ₃ , Sb ₂ O ₃ , others < 0.05%



Heikas Martirosian

Heikas Martirosian, CDT
Technical Director

